

Länge der Achsrohre bezieht sich auf Außenkanten der eingeschweißten Buchsen

Sitzbreite	SPEEDY AF4-4 4all, 4you	SPEEDY 4teen (R50)	SPEEDY AF4-5 4you Ergo, AF5 4all Ergo
SB 26	177mm	187mm	187mm
SB 27	187mm	197mm	197mm
SB 28	197mm	207mm	207mm
SB 29	207mm	217mm	217mm
SB 30	217mm	227mm	227mm
SB 31	227mm	237mm	237mm
SB 32	237mm	247mm	247mm
SB 33	247mm	257mm	257mm
SB 34	257mm	267mm	267mm
SB 35	267mm	277mm	277mm
SB 36	277mm	287mm	287mm
SB 37	287mm	297mm	297mm
SB 38	297mm	307mm	307mm
SB 39	307mm	317mm	317mm
SB 40	317mm	327mm	327mm
SB 41	327mm	337mm	337mm
SB 42	337mm	347mm	347mm
SB 43	347mm	357mm	357mm
SB 44	357mm	367mm	367mm
SB 45	367mm	377mm	377mm
SB 46	377mm	387mm	387mm
SB 47	387mm	397mm	397mm
SB 48	397mm	407mm	407mm
SB 49	407mm	417mm	417mm
SB 50	417mm	427mm	427mm
SB 51	427mm	437mm	437mm
SB 52	437mm	447mm	447mm
SB 53	447mm		
SB 54	457mm		
SB 60	517mm		